



JOUAN EB170 Laboratory Warming Machine (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
French production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 300 W,
Warming machine
JOUAN EB170
has a cast aluminum chamber to ensure an even temperature and the vents effectively remove moisture and fumes,
Specifications:
Infinitely adjustable temperature: 30°C - 60°C,
Chamber capacity: 170 L,
Equipped with additional glass door,
Includes:
Screen trays size: 47 x 47 cm (2 pcs.),
Power cord
Dimensions:
Warmers: 65 x 63.5 x 89 cm,
Chambers: 51 x 51 x 62 cm,
Weight: 70 KG,
It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

