



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



AKRON gynecological chair A9932-8-3FS (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
U.K. production,
Power supply: 230 V,
Emergency power supply: through 2 new built-in 12 V / 1.3 Ah batteries,
Frequency: 50/60 Hz,
New medical approved upholstery in gray (upholstery can be replaced, from our color chart for an additional fee),
AKRON A9932-8-3FS
is used for gynecological and urological procedures,
Specifications:
Electrically adjustable with wired foot control:
Height adjustment in the range: 60 - 97 cm,
Trendelenburg tilt adjustment in range: 17°,
Backrest tilt adjustment within: 70°,
Footrests on both sides with free height adjustment and rotation around its axis,
Paper hanger with width: 47 cm,
Double-sided rails for attaching additional accessories,
Freely attachable leg rest segment,
Mobile chair on 4 wheels (with built-in)
Includes:
Wired foot controller,
Power cord,
Double-sided patient bolsters (2 pcs),
Waste bowl with dimensions: 14.5 x 29 cm,
Slide-in leg rest extension,
Dimensions:
Backrest: 62 x 80 cm,
Seat: 64 x 54 cm,
Leg rest segment: 50 x 50 cm,
Maximum load capacity: 180 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

