



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ENT unit ENTERMED ERGOMED 60 (OLYMPUS OTV-S.C, MAJ-554) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Dutch production,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 2300 VA,
- Bulb: OSRAM HLX 15V / 150W,
- **The ERGOMED 60** unit is designed for ENT examinations. Thanks to its simple design and intuitive control panel, it allows for easy examinations,
- **Specification:**
 - Foiled control panel,
 - Two extendable shelves on guides located under the upper tabletop, possibility of configuring aluminum trays according to your preferences,
 - The drawer is two-story with partitions, the partitions can be configured as desired and there is a stand for viewing glasses,
 - 2 drawers with an opening for used tools and a waste container,
 - Rigid optics heater,
 - Saliva ejector – suction system with a jar for disposable inserts with a capacity of 1 liter,
 - Fiber optic cable,
- **Unity equipped with:**
 - **OLYMPUS OTV-SC Camera Processor with MAJ-554 Camera Head:**
 - Reconditioned (used),
 - Technical condition: very good,
 - Visual condition: very good,
 - Real product photos,
 - Production: Japanese,
 - Power supply: 230 V,
 - Frequency: 50/60 Hz,
 - **Specification:**
 - TV type: PAL,
 - White balance function,
 - Intuitive operation,
 - Keyboard in the control panel in: English,
 - **NEOVO X-15AV Monitor:**
 - Reconditioned (used),
 - Technical condition: very good,
 - Visual condition: very good,
 - Real product photos,
 - Production: PRC,
 - The X-15AV monitor has a screen diagonal of 15 inches, the matrix is additionally protected by an optical filter ensuring excellent color reproduction and good cleaning,
 - **Specification:**
 - Screen diagonal: 15" inches,
 - Screen type: TFT TN,
 - Nominal resolution: 1024 x 768 px.,
 - Contrast: 500:1,
 - Brightness: 400,
 - Spot size [mm]: 0.297,
 - Response time: 12 seconds,
 - Horizontal viewing angle: 140°,
 - Vertical viewing angle: 135°,
- **Included in the set:**
 - OLYMPUS Camera Head Model MAJ-554,
 - OLYMPUS Camera Coupler Model AR-T12E,
- **Dimensions:**
 - Total: 54 x 140 x 95 cm,
- Weight: 230 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for a period of 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
 - **The warranty period does not cover light bulbs and other light sources installed in the device,**
- Financing options (Poland only): Installments, Leasing, Loan,

