



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
  
[kw@praiston.pl](mailto:kw@praiston.pl)



**KARL SOTRZ video-laryngoscopy kit (TELEPACK X 20045020, DUOMAT 20321020, 29003PBK) (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 100 watts,  
KARL STORZ 20045020 TelePack X  
is a portable processor with a built-in light source. It is widely used, from the doctor's office to the operating room, offering a significant advantage in ENT procedures,  
Specifications:  
Monitor:  
18.5-inch touchscreen monitor,  
On-screen keyboard,  
FULL HD,  
Edgeless display - ideal for cleaning,  
Easy to use with intuitive menu navigation and intuitive user interface,  
Light source:  
Integrated LED light source,  
Bright and even illumination,  
Automatic light source control,  
Patient data archiving on 50 GB internal memory, on USB storage media or via network integration function,  
Encrypted patient data transfer,  
Patient data management,  
Selectable menus in ENG, GER, ESP, FR, IT, PL, PT, RUS, JP,  
Set equipped with instrumentation:  
DUOMAT suction and irrigation pump 20321020:  
Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
German manufacture,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 21 VA,  
DUOMAT 20321020  
is a powerful suction and irrigation pump is used in laparoscopic surgeries and can also be used in conjunction with KARL STORZ fiberoptics, Irrigation and suction take place with preset, fixed values. Suction-irrigation tubes with valves can be used to regulate irrigation and suction in laparoscopic applications,  
Specifications:  
Suction power: 3.5 l/min,  
Irrigation power: 3.5 l/min,  
Maximum irrigation pressure: 400 mmHg,  
Maximum suction vacuum: 0.7 bar,  
Instrumentation placed on a mobile video laryngoscope tower KARL STORZ 29003PBK,  
Includes:  
Camera head with fiber optic cable  
KARL STORZ PAL DCI II  
,  
Video-laryngoscope spoon  
KARL STORZ 8401A  
(length: 15.5 cm / width: 2.2 cm),  
Video-laryngoscope spoon  
KARL STORZ 8401B  
(length: 15.5 cm / width: 2.2 cm),  
Wired keyboard  
KARL STORZ 20040240RU  
,  
Suction jar  
KARL STORZ  
,  
Water trap  
ARTEMA 100-000080-00  
(7 pcs),  
BNC-cinch adapter  
KARL STORZ 20040083  
,  
Power cord,  
Overall dimensions: 50 x 55 x 130 cm,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for a period of 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,  
The warranty period does not include bulbs and other light sources installed in the device,  
Possibility of financing (Poland only): Installments, Leasing, Loan,

