



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SCAN MODUL anesthesia assistant (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Dutch production,
Mobile anesthesia assistant
SCAN MODULE
Has top and bottom pull-out drawers along with a retractable side shelf and a waste container mounted on the side rail,
Technical specifications:
Retractable side shelf,
Retractable top drawers for accessories,
Free sliding middle drawers,
Side rails for attaching additional accessories,
Attached waste container on side rails,
Mobile on 4 casters (2 braking),
Dimensions:
Top: 66 x 50 cm,
Retractable side top: 38 x 27 cm,
Top pull-out drawer: 54 x 34 x 5 cm,
3x middle pull-out drawer: 54 x 34 x 11 cm,
Lower pull-out drawer: 54 x 34 x 20 cm,
Total: 73 x 56 x 177 cm,
Weight: 60 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

