



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ATMOS C401 surgical suction machine (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
ATMOS C451
is a versatile surgical suction machine designed mainly for operating rooms and intensive care units,
Specifications:
Front panel pressure gauge,
Maximum flow rate: 40 l/min,
Infinitely adjustable vacuum force up to 91 kPa,
Volume of 48 dB,
Operating mode: continuous operation,
Rail-mounted jar,
Mobile on a mobile base equipped with 4 casters including 2 with a brake,
Included:
Set of germicidal filters,
ATMOS (autoclavable) jar with a capacity of 3 liters (2 pcs.),
Power cord,
Suction tube,
Dimensions: 43 x 38 x 102 cm,
Weight: 14 kg,
It has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

