



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



INNOW S405 treatment and rehabilitation table (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production 2016,
Power supply: 230 V,
Upholstery with medical certificate in cream color (it is possible to replace the upholstery with a new one from our color chart at an additional cost),
Treatment table
INNOW
has a range of height adjustment and backrest tilt which allows free and ergonomic,
Technical specifications:
Manual adjustment:
Height in range: 52 - 75 cm,
Headrest tilt: 26°,
Side free retractable rails,
Bumpers around the table,
Infusion fluid rack with free height adjustment,
Mobile table with built-in brakes,
Includes:
IV fluid rack,
Dimensions:
Headrest: 60 x 44 cm,
Bed: 60 x 144 cm,
Overall length: 190 cm,
Weight: 80 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

