



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **HUNTLEIGH AKRON A3321 medical couch (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
U.K. production,  
Upholstery with medical approval in black (possibility to replace the upholstery with a new one, from our color pattern for an additional fee),  
Couch  
AKRON A3321  
is a two-section couch with adjustable back support,  
Specifications:  
Two-segment couch,  
Structure made of steel profiles covered with powder paint,  
Manual back tilt adjustment in the range: 90° / 0° / -20°,  
Manual height adjustment of the sleep surface in the range: 52 - 100 cm,  
Built-in central brake,  
Paper hanger with a width of 55 cm,  
Dimensions:  
Backrest: 65 x 71 cm,  
Backrest: 65 x 125 cm,  
Overall length: 196 cm,  
Maximum load capacity: 180 KG,  
Weight: 70 KG,  
Has a current inspection and is ready for use,  
Issued Technical Passport (Service Report) valid for a period of 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

