



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SMITH & NEPHEW B65035E endoscope cart (Ex-Display)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production, (2020),
Power supply: 230 V,
Frequency: 50/60 Hz,
SMITH & NEPHEW B650035E
is a stable and sturdy cart suitable for carrying endoscopes and delicate medical equipment,
It has:
Tops with dimensions: 62 x 48 cm,
Retractable lower shelf for accessories with dimensions: 50 x 38 x 8 cm,
Side mount for camera head,
Side mounts for foot manipulator,
Handle for free movement of the cart,
Monitor mounting on 1-segment arm with free height and tilt adjustment,
Electrical installation with the possibility of connecting 8 devices,
Hanger for IV fluids with an adjustable range of 72 cm,
Mobile cart placed on 4 braking wheels,
Can be retrofitted with a monitor at an additional cost,
Includes:
IV stand,
Camera head mount,
Monitor mounting kit,
Dimensions:
Shelves: 62 x 48 cm,
Pull-out drawer: 50 x 38 x 8 cm,
Overall: 80 x 70 x 175 cm,
Maximum load capacity of a single shelf: 50 KG,
Weight: 80 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

