



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



EGERTON SOLID hospital bed (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production,
Mattress with medical certificate in green color,
EGERTON SOLID
is a multifunctional hospital bed designed to meet the needs of the entire medical team,
Due to its functions, it is designed for use not only in general care, but also in post-operative care and in the ICU unit,
It facilitates the use of vital function monitoring devices and patient transport,
Specifications:
Infinitely adjustable manual tilt of the backrest within 70°,
Retractable side rails,
Fixed recliner height at 61 cm,
Drip stand height adjustment within: 68 cm,
Mobile bed on 4 braking casters,
Includes:
Drip stand,
Dimensions:
Mattress: 90 x 196 cm,
Overall: 100 x 216 x 80 cm,
Maximum load capacity: 170 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

