



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ANETIC AID QA3 transport bed (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: good,
Real photos of the product,
U.K. production, (2019),
Anetic Aid QA3
is a functional transport bed that can also serve as an operating table thanks to its numerous functions,
Extremely maneuverable, compact, lightweight and durable for maximum patient comfort,
Medical approved mattress in black,
Bed made of double laminated reinforced HPL board,
Has visible signs of use,
Adjustable in many planes (hydraulically and through gas springs):
Bed height in the range: 63.5 - 93 cm,
Trendelenburg position in the range: 14°,
Backrest angle by means of in the range of 0° to 90°,
Features:
X-RAY translucency,
Folding handrails,
Retractable IV stands in the handrail, Height adjustable from: 140 - 180 cm,
Plastic base, with space for a cylinder,
Round rubber bumpers on the corners,
Mobile on four casters with diameter: 15 cm,
Central brake and the possibility of directional travel with the help of the extendable fifth wheel,
Dimensions:
Mattress: 70 x 202 x 10 cm,
General: 76 x 216 x 84.5/114 cm,
Maximum load capacity: 250 kg,
Weight: 135 KG,
Cat 03,
Has a current Technical Passport issued,
Warranty:
6 months for domestic market - Poland,
3 months for the international market,
Possibility of financing (Poland only): Installments, Leasing, Loan,

