



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **ESSILOR DIGITAL CRP digital pupilometer (Reconditioned)**

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

P.R.C. production,

Power supply: 2 x AA battery (1.5 V),

CRP Digital Pupilometer:

Corneal Reflectance Pupilometer is used for precise measurement of pupil distance and its numerical recording,

The Pupilometer measures monocular and binocular pupil distance for all distances from 35 cm to infinity,

Monocular pupil differences can be measured by closing the eye if the patient's binocular vision is abnormal,

The pupil differences between monocular and total pupils (48 to 77 mm) are measured to the nearest 0.5 mm,

Features automatic illumination and digital display,

Dimensions: 25 x 16 x 6 cm,

Weight: 0.7 KG,

It has a current inspection and is ready to work,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

