



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS UES-40 SURGMASTER surgical diathermy (Reconditioned)

Product temporarily unavailable
Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Japanese production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 1380 W,
Technical specifications:
Monopolar - max. 300 W/300Ω,
Cutting - URO / PURE / BLEND,
Coagulation - COAG.1 / COAG.2 / SPRAY,
Maximum power:
PURE: 300 W,
BLEND: 250 W,
URO: 300 W,
COAG.1: 120 W,
COAG.2: 120 W,
SPRAY: 120 W
Bipolar, 90 W/100 Ω max,
Cutting - PURE,
Coagulation - SOFT1, SOFT2, HARD,
Maximum power:
PURE: 90 W,
SOFT.1: 90 W,
SOFT.2: 90 W,
HARD:
L1: 80 W,
L2: 120 W,
L3: 160 W,
Saline, max. 320 W/100 Ω
(resection in saline, e.g., during TUR),
Cutting - PURE, BLEND,
Coagulation - COAG.1, COAG.2,
Maximum power:
PURE: 320 W,
BLEND: 320 W,
COAG.1: 200 W,
COAG.2: 80 W,
Keypad on control panel in language: English,
Included:
User manual in English (PDF),
Power cord,
Factory new cable with ERBE clip for neutral electrodes
- No electrodes included, can be retrofitted at additional cost,
Factory new monopolar electrode cable with ERGO PEN BOWA JackKNIFE holder - 2 buttons, 2.4 mm shaft (No electrodes included, can be retrofitted at additional cost),
Factory new bipolar cable, BOWA forceps
(can be retrofitted with bipolar tweezers at additional cost),
Wired foot attachment,
Dimensions: 41 x 35 x 15.5 cm,
Weight: 15 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

