



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GRAND SEICO GL-7000B Dioptrimeter (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Japanese production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 40 VA,
Specifications:
Built-in 5.7" inch color LCD touchscreen display,
Measuring range - Sphere (D): -25 - 25 (gradual adjustment: 0.01/0.12/0.25),
Measuring range - Cylinder (D): 0 - 10 (gradual adjustment: 0.01/0.12/0.25),
Measuring range - Prism (D): 0 - 10 (gradual adjustment: 0.01/0.12/0.25),
Wavelength: 660 nm,
Automatic recognition of progressive and multifocal lenses,
Prismatic test lenses,
UV transmittance measurement,
Type of test lenses:
Glass lenses,
Hard / soft contact lenses,
Maximum lens size:
100 mm,
Selectable menu in language: ENG, CHIN,
Included:
English user manual (PDF),
Power cord,
Dimensions: 20 x 39.5 x 36.5 cm,
Weight: 6 KG,
Has an up-to-date inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

