



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CONMED 10K irrigation pump (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Manufactured by U.S.A., (2015),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 65 VA,
Technical specifications:
Maximum pressure: 150 mmHg,
Maximum flow rate: 2000 ml/min,
Maintains stretch even in the most demanding situations,
Intuitive one-handed cartridge loading,
Includes:
Power cord,
Dimensions: 25 x 20 x 18 cm,
Weight: 3 KG.
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

