



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



TOMEY TL-2000A Lens Meter (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Japan,
Power supply: 230 V,
Frequency: 50/60 Hz,
Features:
Very easy to use,
Automatic measurement and storage of measurement results,
Automatic detection of progressive lenses,
Contact lens measurement,
Built-in thermal printer: (paper width 38 mm),
Built-in display: Backlit, monochrome, liquid crystal,
Features:
Measuring range:
Sphere: from -25 to +25D,
Cylinder: from -9.99 to +9.99D,
Axis: from 0° to 180°,
Addition: from 0 to 9.99D,
Prism: from 0 to 10D,
Measurement step:
every 0.01/0.12/0.25D,
Measurement modes:
Cylinder: +/-
, mixed,
Prism:
rectangular, polar,
Lens supported:
Standard, bifocals, trifocals, progressive, high-index, contact lenses,
Lens supported diameters:
from 20 to 100 mm,
ABBE numbers:
30 / 35 / 40 / 45 / 50 / 55/ 60,
Menu language: English,
Included:
Power cord,
Dimensions: 20 x 29 x 42 cm,
Weight: 10 KG,
Updated inspection and ready to use,
Issued Technical Passport (Service Report) valid for 12-month warranty
Warranty:
6 months for the domestic market (Poland),
3 months for the international market
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets
Financing options (Poland only): Installments, Leasing, Loan

