



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



KARL STORZ 11301 BN1 fiberscope (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany,
Specifications:
Viewing angle: 120°
Upward tilt: 140°
Downward tilt: 140°
Distal outer diameter: 5.2 mm,
Working length: 65 cm,
Working channel diameter: 2.3 mm,
Included:
Dedicated KARL STORZ H2O2 transport case,
Case key (2 pcs),
Weight: 0.35 kg,
Updated inspection and ready to use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee up to 12 months for domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

