



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



POLLOGEN GENE0 skin care device (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Produced in Israel (2017),
Power supply: 230 V,
Frequency: 50/60 Hz,
POLLOGEN GeneO+

- is a Super-Facial platform that represents the latest technology in facial skin care. It utilizes "OxyGeneo" technology, providing excellent skin nourishment, oxygenation, and long-lasting skin firming, as well as proven anti-aging results.

Features:

OxyGeneo offers a unique approach to facial treatments. It utilizes natural physiological processes to increase oxygen levels in the skin and achieve exceptional results. Unlike other technologies, OxyGeneo utilizes natural oxygen from within the skin, ensuring effective treatment and immediate results.

The resulting CO₂-rich environment initiates a physiological response known as the "Bohr effect," leading to increased tissue oxygenation, resulting in:

Increased blood flow,
Increased oxygen release to the area,
Increased skin cell metabolism.

Specifications:

10.2" color LCD touchscreen display,
Side head holder,
Menu language selection: GER, ENG, ESP, FR, IT, NED, PL, RUM, SK, TUR, CS, RUS, CHIN,

Includes:

English user manual (PDF),
Power cord,
GENOE0 25002020 head

Dimensions: 38 x 28 x 38 cm,

Weight: 4.4 kg KG,

It has a current inspection and is ready for operation.

A Technical Passport (Service Report) is issued, valid for 12 months.

Warranty:

6 months for the domestic market (Poland),

3 months for the international market.

The warranty can be extended for an additional fee to 12 months for the domestic and international markets.

Financing options (Poland only): Installments, Leasing, Loan.

