



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SML 32 / 250 ZALMED laboratory dryer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Poland,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 1000 W
Capacity: 32 l,
Enables drying, heating, and sterilization with dry, hot air from ambient temperature to 250°C in an air atmosphere,
Perfectly suited for equipping pharmacies as a dryer and as an air sterilizer (PZH certified) used for drying and sterilizing glass, porcelain, ceramic, and metal products,
Specifications:
Temperature control: From 0°C to 250°C,
9-level drying power adjustment,
Includes:
Shelf dimensions: 23 x 40 cm (1 pc.),
Shelf dimensions: 27 x 40 cm (1 pc.) pcs.),
Overall dimensions: 50 x 38 x 60 cm,
Weight: 32 kg,
Updated inspection and ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

