



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Video - EVIS EXERA II GIF gastroscope - Q165 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Japan,
Power supply: 230 V,
Frequency: 50/60 Hz,
The Olympus GIF-Q165
is a video gastroscope that produces high-quality images.
The outer diameter of the endoscopic probe tip and the outer diameter of the endoscopic probe tube provide excellent insertion options.
The Scope ID function stores information in a built-in memory system and displays it on the monitor.
The wide field of view allows for precise observation of the selected area.
The 2.8 mm diameter endoscope channel allows for the use of a wide range of therapeutic accessories.
Convenient button access and user-programmable switches simplify operation.
The 4-way bendable endoscope tip allows for a complete and comprehensive examination of the upper gastrointestinal tract.
Technical parameters:
Optical system:
Field of view: 140°.
Direction Viewing Distance: Front View,
Depth of Field: 3 - 100 mm,
Endoscopic Probe:
Outer Diameter of the Endoscopic Probe Tip: 9.2 mm,
Outer Diameter of the Endoscopic Probe Tube: 9.2 mm,
Working Length: 1030 mm,
Biopsy Channel:
Inner Diameter of the Channel: 2.8 mm,
Minimum Viewing Distance: 3.0 mm from the Endoscopic Probe Tip,
Bending Tip:
210° Up,
90° Down,
100° Right,
100° Left,
Total Length: 1345 mm,
Fully Compatible: CV-160 / CV-145 / CV-140,
Included:
OLYMPUS MH-443 Suction Valve,
OLYMPUS MH-443 Air-Water Valve 438,
Cover,
OLYMPUS VB-301 Biopsy Valve,
Biopsy Valve,
Dedicated transport case with dimensions: 59.5 x 43.5 x 13 cm,
Cat 09,
Updated inspection and ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan

