



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BLOODY 04 blood collection chair (New)

The new, Polish-made BLOODY 04 blood collection chair was designed with patient comfort and the ergonomics of medical staff in mind during procedures such as blood collection, IV drip administration, and injections. The chair is intended for use in medical facilities such as laboratories, clinics, hospitals, and diagnostic centers.

It allows for comfortable and safe patient positioning during blood draws or other minor outpatient procedures.

The stable structure made of powder-coated steel ensures durability and ease of disinfection.

The seat and footrest are adjustable via a gas spring to the Terndelenburg position, which helps in proper patient positioning in the event of fainting.

Upholstery in your choice of color is available free of charge, except for the "PRESTIGE" colors.

The chair can be equipped with:

Standard headrest,

Footrest cover,

Metal instrument tray,

Technical Specifications:

Height 100 cm,

Width 50 cm, Total 71 cm,

Depth 81 cm,

Total length when extended 143 cm,

Seat height 50 cm cm,

Own weight: 30 kg,

Maximum load capacity: 150 kg,

Headrest height adjustment: 15 cm,

Armrest height adjustment: 25 cm,

Backrest angle adjustment: 0° to 85°,

Footrest angle adjustment: 5° to -90°,

Valid Technical Passport issued,

24-month warranty on the structure,

