



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



---

### **CAZANDER MEDICAL 361 anesthesia assistant (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
CAZANDER MEDICAL 361 ASSISTANT  
is stable and sturdy and also suitable for anesthesia procedures,  
Specifications:  
Tilting top drawers for accessories (8 pcs),  
Upper pull-out drawers (2 pcs),  
Sliding lower drawers for accessories (4 pcs),  
Side rail suitable for attaching additional accessories,  
Attached stainless steel upper shelves (4 pcs),  
Attached side basket,  
Attached side handle for transporting the assistant,  
Cart placed on 4 casters (2 braking),  
Dimensions:  
Upper top: 80 x 40 cm,  
Stainless steel upper shelves: 14 x 58 cm (4 pcs),  
Top drawers: 70 x 38 x 10 cm (2 pcs),  
Bottom drawers: 30 x 40 x 16 cm (4 pcs.),  
Weight: 100 KG,  
Has an up-to-date inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

