



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Gynecological and treatment chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Power supply: 230 V,
Frequency: 50/60 Hz,
Medically certified blue upholstery (upholstery can be replaced with a new one from our color chart for an additional fee),
Technical specifications:
Electric height adjustment: 48 - 78 cm,
Manual backrest adjustment: 90°,
Manual legrest adjustment: 90°,
Sliding waste pan: 17 x 26 cm,
Bilateral knee rests with free manual height adjustment,
Free manual headrest height adjustment,
Included:
Wired remote control,
Dimensions:
Headrest: 32 x 16 cm,
Backrest: 60 x 65 cm cm,
Seat: 60 x 46 cm,
Legrest: 60 x 46 cm,
Weight: 51 kg,
Updated inspection and ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international market,
Financing options (Poland only): Installments, Leasing, Loan,

