



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical chair (for patient transport) HILL-ROM ANATOMIE SM647B (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in France,
Navy blue upholstery with medical certificate – visible signs of use with minor damage to the seat (upholstery can be replaced with a new one from our color chart for an additional fee),
Specifications:
Hydraulic height adjustment in the range of: 470 – 700 mm,
Smooth backrest adjustment in the range of 24° (gas spring),
Smooth seat angle adjustment (so-called Cradle),
Adjustable, retractable and extendable armrests with a range of 200 mm,
Headrest adjustable in 3 positions,
Retractable and retractable footrest,
Three-position brake - wheel lock control,
Movement handles,
Washable upholstery,
Infusion stand height: 1500 mm,
Dimensions:
Overall: 680 x 1300 x 1330 mm,
Backrest: 460 x 550 mm,
Seat: 470 x 450 mm,
Maximum weight capacity: 135 kg,
Weight: 58 kg,
Up-to-date inspection and ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty: 6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

