



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDIFA doctor's stool (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany,
Medically certified black upholstery (upholstery can be replaced with a new one from our color chart for an additional fee),
Height adjustable with a gas spring,
Infinitely adjustable seat height from 42 to 55 cm,
Base diameter: 60 cm,
Mobile on double wheels,
Seat can rotate 360°,
Dimensions:
Seat diameter \varnothing : 34 cm,
Weight: 10 kg,
Updated inspection and ready to use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market international,
Possibility of extending the warranty for up to 12 months for domestic and international markets for an additional fee,
Financing options (Poland only): Installments, Leasing, Loan,

