



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Doctor's chair (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- **Upholstery with medical certificate in black (possibility of replacing the upholstery with a new one from our color sample for an additional fee),**
- Height adjustable by means of a gas cylinder,
- Smooth seat height adjustment in the range: 50 – 68 cm,
- Base diameter: 50 cm,
- Mobile on double wheels,
- The seat can be rotated around its own axis in the range of 360°,
- Dimensions:
 - Seat diameter ø: 32 cm,
- Weight: 10 KG,
- Cat 02,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

