



HETTICH ROTINA 46R laboratory centrifuge (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Made in Germany,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 1350 W,

The HETTICH ROTINA 46R centrifuge is a refrigerated tabletop centrifuge with a maximum capacity of 4 x 400 ml. This centrifuge is designed to separate substances or mixtures with a density of up to 1.2 kg/dm³.

Specifications:

Maximum rotor speed: 4000 rpm,

Maximum capacity: 4x 400 ml,

Maximum density: 1.2 kg/dm³,

Maximum speed: 11,000 rpm,

Relative centrifugal force: 13,800,

Noise level: ≤ 66 dB(A),

Temperature range: -5°C - +40°C,

Includes:

Power cord,

HETTICH ROTINA 4624 rotor,

4x HETTICH ROTINA 4630 jugs,

2x HETTICH ROTINA 4639 adapters,

2x HETTICH ROTINA 4638 adapters,

Overall dimensions: 81.5 x 69 x 41 cm,

Weight: 88 kg,

It has a current inspection and is ready to work,

A Technical Passport (Service Report) is issued, valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility to extend the warranty for an additional fee to 12 months for the domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan,

