



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Gynecological chair MAQUET RADIUS 1557.00AJP2 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 290 VA,
Maquet Radius
is used for gynecological and urological examinations and procedures,
Medically certified blue upholstery (upholstery can be replaced with a new one from our color chart for an additional fee),
Electromechanical adjustment via handheld remote control and foot controller:
Bed height range: 56 - 105 cm,
Tilt angle:
Back section range: (- 11°) / 0° / 57°,
Seat angle range: 25°,
When adjusting the back section angle, the tilt angle automatically adjusts seat,
Storage compartment for a 38 cm wide roll of backing paper,
Cut-out in the back section upholstery for backing paper,
3 pre-programmed buttons on the handheld remote,
STOP button on the handheld remote to stop the programmed button,
Electric recline of the seat in a lying position,
Free adjustment of the footrest by 30° horizontally,
Includes:
Instruction manual in Polish (PDF),
Wired foot controller,
Wired handheld remote,
Power cord,
Waste bowl measuring 28 x 18 cm,
Headrest,
Dimensions:
Back section: 91 x 61 cm,
Seat section: 61 x 31 cm,
Length in the unfolded position: 175 cm,
Overall dimensions: 155 x 76 x 138/168 cm,
Maximum Load capacity: 180 kg,
Weight: 94 kg,
It has a current inspection and is ready for operation,
A Technical Passport (Service Report) is issued, valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

