



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



WS-TECH LM-01 medical couch (Ex-Display)

New (ex-display),
Actual product photos,
Made in Poland (2024),
New PRESTIGE upholstery with medical certification in light brown (upholstery can be replaced from our color chart for an additional fee),
The LM-01 stationary couch is a classic model designed for performing physical examinations, administering injections, or removing stitches.
This stationary couch is perfect for medical and physiotherapy offices, as well as professional massage and tattoo parlors.
It features a 45° adjustable headrest height.
Built-in paper holder.
Dimensions:
Headrest: 70 x 40 cm,
Total: 70 x 190 x 70 cm,
Maximum load capacity: 150 kg,
Weight: 34 kg.
It has a current inspection and is ready for use.
A Technical Passport (Service Report) is issued and is valid for 12 months.
Warranty:
6 months for the domestic market (Poland),
3 months for the international market.
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international market.
Financing options (Poland only): Installments, Leasing, Loan.

