



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PHILIPS HD11XE ultrasound machine (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Made in the U.S.A.,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 1000 VA,

The

PHILIPS HD11XE ultrasound machine

is an ultrasound machine tailored to the needs of every user and provides a complete solution for high-resolution digital imaging.

PHILIPS has equipped the machine with a broadband beamforming system and automatic image optimization tools. This solution makes ultrasound ideal for public and private hospitals, clinics, and specialized practices.

Display: 17-inch LCD (rotatable and adjustable).

Built-in thermal printer

MITSUBISHI P95

Print formats: 100 x 75 mm / 100 x 130 mm.

Imaging technologies:

SonoCT:

Multiple-line-of-sight image fusion technology (up to 9) for improved clarity.

XRES:

Adaptive image processing reduces artifacts and improves tissue boundary definition.

Harmonic Imaging:

Harmonic imaging enhances contrast.

Doppler modes:

Color Doppler, Power Doppler, Pulsed Wave (PW), Continuous Wave (CW) with Adaptive Doppler technology (weak signal enhancement, noise reduction).

Ergonomics and portability:

Adjustable panel Control panel and monitor, lockable wheels, easy to transport,

Data management:

Basic archiving functions on CD/DVD, creating reports with embedded images, on-screen image thumbnails,

Clinical applications:

Basic: Abdominal, obstetrics/gynecology (OB/GYN), small organs, vascular,

Menu language: English,

Included:

PHILIPS C5-2 convex probe,

PHILIPS L12-3 linear probe,

PHILIPS S4-2 sector probe,

PHILIPS C8-4V endovaginal probe,

Instruction manual in English (PDF),

Power cord,

The device has visible signs of use that do not affect comfort of use,

Dimensions: 53 x 110 x 130 / 148 cm,

Weight: 100 kg,

Valid inspection and ready for use,

Technical Passport issued (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Warranty can be extended for an additional fee to 12 months for the domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan,

