



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MIDMARK PROMOTAL 3181 treatment and cosmetic chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
French production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 230 W,
Upholstery with medical certificate in gray color (it is possible to replace the upholstery with a new one from our color pattern for an additional fee),

Treatment and cosmetic chair

MIDMARK PROMOTAL 3181

provides free adjustment with positioning in a wide range,

The chair has excellent mobility due to which it is adapted for free movement around medical institutions,

Technical specifications:

Adjustment by wired handheld remote control:

Seat height adjustment in the range: 56 - 80 cm,

Backrest adjustment in the range: -7° / 0° / 74°,

Adjustment of the leg rest in the range: -57° / 0° / 5°,

Tilt to Trendelenburg position,

The remote control has a pre-programmed button with the setting of the chair to the initial position,

Mobile chair on 4 wheels with built-in brake,

Freely adjustable patient armrests on both sides,

Freely adjustable drip rack with both sides,

Includes:

Patient headrest,

2 pcs. - armrest,

IV hanger

,
Dimensions:

Backrests: 58 x 80 cm,

Seats: 58 x 53 cm,

Leg rests: 58 x 46 cm,

Overall: 95 x 115 x 145/169 cm,

Maximum load capacity: 135 KG,

Weight: 95 KG

It has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

