



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ULTRAVIOL NBV30PLW germicidal lamp (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Made in Poland,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 33 W,

The device is designed to improve and maintain the microbiological cleanliness of rooms. It emits UV-C radiation.

Irreversibly deactivates bacteria, viruses, molds, fungi, and all other microorganisms. It has a reflector made of very high-quality reflective aluminum.

Specifications:

UV-C radiation intensity at a distance of 1 m: 2.3 W/m²

Lamp life: 9000 h

Disinfected area: 12-15 m²

Lamp rotation (adjustable beam angle): 270°

Dimensions:

Dome: 96 x 8.5 x 13.5 cm

Arm length: 180 cm

Weight:

Dome weight: 3 kg

It has a current inspection and is ready for operation

A Technical Passport (Service Report) is issued, valid for 12 months

Warranty:

6 months for the domestic market (Poland),

3 months for the international market

Possibility to extend the warranty for an additional fee 12 months for domestic and international markets.

The warranty period does not cover light bulbs and other light sources installed in the device.

Financing options (Poland only): Installments, Leasing, Loan.

