



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ARNOLD CAE Mobile Chair Scale (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in France,
Power supply: 230 V,
Backup power supply: via 1 new built-in 6 V / 1.2 Ah battery,
Frequency: 50/60 Hz,
ARNOLD CAE Mobile Chair Scale
characterized by ease of use and comfort for the patient, with 2 built-in wheels for easy transport,
Specifications:
Retractable armrests on both sides,
Retractable backrest,
Extendable legrest with a range of 17 cm,
Built-in electric patient scale:
Weighing accuracy: 100 g,
Weighing position: standing / sitting,
Digital LCD display,
Dimensions:
Backrest: 44 x 52 cm,
Total: 59 x 59 x 105 cm,
Maximum load capacity: 200 kg,
Weight: 30 kg,
It has a current inspection and is ready for operation,
A Technical Passport (Service Report) is issued, valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

