



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MDT MD-2N Ophthalmic Table (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Poland,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 80 VA,
2 built-in 230 V sockets for powering additional accessories,
Electric height adjustment in the range: 74 - 94 cm,
Included:
Power cord,
Dimensions:
Tabletop: 42 x 95 cm,
Weight: 35 kg,
Valid inspection and ready to use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic market and international,
Financing options (Poland only): Installments, Leasing, Loan,

