



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



TOPCON CT-80A non-contact tonometer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Japan,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 80 VA,
Topcon CT - 80A
is a non-contact tonometer with a reduced air blast force, which reduces patient discomfort during the examination.
Non-contact,
Measuring range: 0 - 60 mmHg,
Software: 2.1.,
Features:
5" monitor,
Device menu in English,
Triple safety function:
Adjustable safety stopper,
"Too close" message,
Audible signal warning when the device is too close to the cornea,
Measurement mode: Automatic start with automatic 3D measurement alignment or manual (up to selection),
Automatic pupil targeting,
Electrical bowl adjustment,
Manually adjustable chin rest,
Built-in RS-232 port,
Thermal printer with a paper width of 5.7 cm,
Menu language selection: English,
Included:
Power cord,
Dimensions:
27.5 x 46 x 47 cm
Weight: 19 kg,
Updated inspection and ready to use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international market,
Financing options (Poland only): Installments, Leasing, Loan,

