



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ANETIC AID QA3 transport bed (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- UK Production, (2015),
- Anetic Aid QA3 is a functional transport bed that can also be used as an operating table thanks to its numerous functions,
- Extremely agile, compact, lightweight and durable, ensuring maximum patient comfort,
- Medically certified mattress in black,
- Bed base made of double laminated reinforced HPL board,
- Has visible signs of use,
- **Multi-level adjustment (hydraulic and via gas springs):**
 - Bed height range: 63.5 - 93 cm,
 - Trendelenburg position in the range: 14°,
 - Backrest angle adjustable from 0° to 90°,
- **It has:**
- X-RAY transparency through a full-size X-ray platform,
 - Folding handrails,
 - IV poles hidden in the railing, adjustable in height from 140 to 180 cm,
 - Plastic base with space for a bottle,
 - Round rubber bumpers on the corners,
- Mobile on four wheels with a diameter of: 15 cm,
- Central brake and the possibility of directional driving using a retractable fifth wheel,
- Dimensions:
 - Mattress: 70 x 202 x 10 cm,
 - Overall: 76 x 216 x 84.5/114 cm,
- Maximum load capacity: 250 KG,
- Weight: 135 KG,
- Cat 04,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

