



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **STOLTER MOBILO transport trolley (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in Poland (2014),  
2-section,  
Waterproof mattress with medical certificate, beige,  
Specifications:  
Manual, gradual headrest adjustment,  
Manual bed height adjustment using a crank within the range: 65 – 87 cm,  
Double-sided bumpers,  
50 cm wide paper holder,  
Transparent tabletop suitable for X-rays,  
Double-sided retractable rails,  
Mobile trolley mounted on 4 braking wheels,  
Includes:  
Infusion rack,  
Mattress dimensions: 62 x 182 cm,  
Weight: 100 kg,  
Updated inspection and ready for use,  
Issued Technical Passport (Service Report) valid for 12-month warranty  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market  
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets  
Financing options (Poland only): Installments, Leasing, Loan

