



MAQUET 1002.72A0-04 headrest with MAQUET 1002.64A0 articulated handle (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: good,
- **Real product photos,**
- German production,
- **MAQUET 1002.72A0-04** is a two-segment, divided headrest,
- **Specification:**
 - It is used to support the patient's head in a lying position, both on the back and on the stomach,
 - Manual adjustment of the headrest in various positions,
 - Headrest width adjustment (inside) in the range: 9 - 18 cm, each side separately,
 - 2 sections of the headrest are adjustable independently of each other,
 - Suitable for use with **MAQUET** surgical table : 1130 / 1131 / 1132 / 1140 / 1150,
- Overall dimensions: 45 x 25 x 34 cm
- Weight: 7 KG,
- Cat 02,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

