



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CONVEX GE RAB4-8L head (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Manufactured by P.R.C.,
GE RAB4-8L

is a convex probe with a frequency range of 4.0 – 8.0 MHz. It is a 4-dimensional, convex volumetric transducer used in obstetrics and gynecology, abdominal, and urological ultrasound.

Frequency range: 4.0 – 8.0 MHz,

It has a current inspection and is ready to operate.

A Technical Passport (Service Report) valid for 12 months is issued.

Warranty:

6 months for the domestic market (Poland),

3 months for the international market.

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets.

Financing options (Poland only): Installments, Leasing, Loan.

