



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SIEMENS EC9-4W endovaginal head (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Manufactured in the U.S.A. (2015),
Siemens
EC9-4W
is an intra-abdominal ultrasound transducer operating in the frequency range from 4.0 to 9.0 MHz,
Specifications:
Frequency range: 4.0 – 9.0 MHz,
Array type: intracavitary,
Imaging modes: 2D, Color Doppler, Power Doppler, Harmonic imaging,
Field of view: 170°,
Cable length: 2.1 meters,
Features:
Frequency range:
from 4.0 to 9.0 MHz, suitable for detailed imaging of pelvic structures,
Intracavitary design:
allows you to get closer to the area of interest, increasing image resolution,
Ergonomic shape:
Lightweight and User-friendly design reduces operator fatigue.
Durable construction:
Designed to withstand routine clinical use and disinfection processes.
Applications:
Obstetric imaging:
Detailed visualization of fetal development and maternal health.
Gynecological imaging:
Assessment of pelvic organs and reproductive health.
Urological imaging:
Assessment of urological structures and conditions.
Fertility assessment:
Assistance in fertility assessment and procedures.
Updated inspection and ready for operation.
Issued Technical Passport (Service Report) valid for 12 months.
Warranty:
6 months for the domestic market (Poland),
3 months for the international market.
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international markets.
Financing options (Poland only): Installments, Leasing, Loan.

