



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ANETIC AID QA4 transport cart (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- UK production,
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in batteries 12 V / 4.5 Ah,
- Frequency: 50/60 Hz,
- **Anetic Aid QA4** is a functional transport trolley that can be used as an operating table thanks to its numerous functions,
- **Upholstery with medical certificate in black (possibility of replacing the upholstery with a new one for an additional fee),**
- Bed base made of double laminated reinforced HPL board,
- The mattress is covered with a waterproof, antibacterial coating,
- **Manual adjustment:**
 - Headrest in the range of (-30)° / 0° / 25°,
 - Lateral tilt in the range of 12° using a crank,
- **Electric control via remote control:**
 - Heights ranging from 68 to 108 cm,
 - Backrest inclination angle within 65°,
 - Longitudinal travel in the range of 45 cm,
 - Changes in the patient's lying direction relative to the table column: Trendelenburg position within 20° and anti-Trendelenburg position within 12°,
- **Equipped with:**
 - Handles on both sides of the headboard for convenient bed guidance,
 - Rails for mounting additional equipment,
 - Groove in the base for an oxygen cylinder,
 - 2 convenient brakes on both sides of the bed,
 - Directional wheel with lever,
 - 2 fender discs,
 - Emergency handle allowing the user to immediately lower the backrest in case of an emergency,
 - Retractable barriers on both sides,
 - Foldable IV stand,
- Table adapted to C-arm,
- Possibility of disassembling the leg section and headrest,
- The table has X-RAY transparency,
- Mobile, on 4 double wheels with a diameter of 15 cm,
- **Included in the set:**
 - User manual in English (PDF),
 - Power cord,
 - Wired remote control,
- Dimensions
 - Overall: 211 x 80 x 68/108 cm,
 - Bed base: 200 x 59 x 8.5 cm,
- Maximum load capacity:
 - Beds: 190 kg,
 - Head segment: 50 kg,
 - Leg segment: 25 kg,
- Weight: 153 KG,
- Cat 04,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

