



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **FAMED ŻYWIEC LP-01.5 rehabilitation bed (Ex-Display)**

New (ex-display),

Actual product photos,

Made in Poland,

The FAMED ŻYWIEC LP-01.5 is a comfortable rehabilitation bed with manual adjustment.

Specifications:

Fixed bed height: 50 cm,

Gradual headrest adjustment in the range of: 0° / 68°,

Bump stop on each side of the bed,

Solid and durable construction,

Mobile bed placed on 4 wheels (2 with brakes),

Overall dimensions: 80 x 180 x 95 cm,

Weight: 25 kg,

Updated inspection and ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility to extend the warranty for an additional fee to 12 months for the domestic market and international,

Financing options (Poland only): Installments, Leasing, Loan,

