



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FAMED ŻYWIEC LR-07.0 rehabilitation bed (Ex-Display)

New (ex-display),

Actual product photos,

Made in Poland,

The FAMED ŻYWIEC LR-07.0

is a comfortable rehabilitation bed with manual adjustment.

Technical specifications:

Manual adjustment:

Height: 55 - 88 cm,

Backrest adjustment: 0° / 67°,

Legrest adjustment: 0° / - 40°,

Trendelenburg adjustment: 25°,

Reverse Trendelenburg adjustment: 12°,

Body bumpers on each side of the bed,

By-sided foldable side rails,

Built-in central brake and guidance system,

Dimensions:

Mattress: 65 x 192 cm,

Weight: 70 kg,

Up-to-date inspection and ready for use.

Issued Technical Passport (Service Report) valid for a period of 12-month warranty

Warranty:

6 months for the domestic market (Poland),

3 months for the international market

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets

Financing options (Poland only): Installments, Leasing, Loan

