



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **CONVEX GE AB2-7 ultrasound head (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Manufactured by P.R.C.,  
GE AB2-7

is a convex probe with a frequency range of 2.0 - 7.0 MHz. It supports the following applications: pediatrics, obstetrics, gynecology, urology, and the musculoskeletal system.

Frequency range: 2.0 - 7.0 MHz,

It has a current inspection and is ready to operate.

A Technical Passport (Service Report) valid for 12 months is issued.

Warranty:

6 months for the domestic market (Poland),

3 months for the international market.

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets.

Financing options (Poland only): Installments, Leasing, Loan.

