



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **STOLTER USB-04 diagnostic and treatment table (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in Poland,  
Medically certified yellow upholstery (upholstery can be replaced with a new one from our color chart for an additional fee),  
USB-04

is a universal examination table for medical, treatment, gynecological, and urological healthcare facilities.

The table is made of steel profiles with a powder coating that is resistant to mechanical damage.

It has four seamlessly upholstered bed segments.

The table has manual adjustment:

Headrest tilt angle: 0 - 45°,

Backrest tilt angle: 0 - 75°,

Seat tilt angle: 0 - 15°,

Trendelenburg position: 0 - 10° and reverse Trendelenburg position: 0 - 15° with gas spring support,

Base height adjustment via foot lever: 75 - 104 cm,

Table equipped with:

Stainless steel rails on both sides for mounting additional accessories,

Stainless steel pull-out tray with dimensions: 30 x 24 x 6.5 cm,

Paper holder with width: 51 cm,

2 patient footrests,

Dimensions:

Headrest: 38 x 32 cm,

2-segment backrest: 62 x 162 cm,

Seat: 62 x 36 cm,

Total length: 196 cm,

Weight: 80 kg,

Up-to-date inspection and ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty: 6 months for the domestic market (Poland),

3 months for the international market,

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan,

