



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STOLTER USB-04 gynecological chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Poland,
Medically certified yellow upholstery (upholstery can be replaced with a new one from our color chart for an additional fee),
USB-04 is a universal examination table for medical, treatment, gynecological, and urological healthcare facilities.
The table is made of steel profiles with a powder coating resistant to mechanical damage.
It has four seamlessly upholstered bed segments.
The table has manual adjustment:
Headrest tilt angle: 0 - 45°,
Backrest tilt angle: 0 - 75°,
Seat tilt angle: 0 - 15°,
Trendelenburg position: 0 - 10° and reverse Trendelenburg position: 0 - 15° with gas spring support,
Base height adjustment via foot lever: 75 - 104 cm,
Table equipped with:
Stainless steel rails on both sides for mounting additional accessories,
Stainless steel pull-out tray with dimensions: 30 x 24 x 6.5 cm,
Paper holder with width: 51 cm,
2 patient footrests,
Dimensions:
Headrest: 38 x 32 cm,
2-segment backrest: 62 x 162 cm,
Seat: 62 x 36 cm,
Total length: 196 cm,
Weight: 80 kg,
Up-to-date inspection and ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty: 6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

