



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ENRAF NONIUS MANULEKT-U treatment and diagnostic table (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in the Netherlands,
Power supply: 230 V,
Frequency: 50/60 Hz,
Medically certified black upholstery (upholstery can be replaced with a new one from our color chart for an additional fee),
Electrically adjustable table height using a hydraulic cylinder,
Manual position adjustment,
Retractable wheels and adjustable feet provide excellent stability (RWD),
Frame construction designed for easy use as a mobile lift,
Powder-coated frame with an antibacterial coating,
Technical specifications:
Electric height adjustment in the range: 46 - 104 cm,
Manual adjustment:
Headrest adjustment: +33° / 0° / -20°,
Backrest adjustment: +80° / 0°,
Legrest adjustment: +24° / 0°,
Reflective backrest/seat and seat/legrest positions
REFLEX
Dimensions:
Headrest: 67 x 35 cm,
Backrest: 60 x 35 cm,
Seat: 48 x 35 cm,
Legrest: 60 x 35 cm,
Total length: 205 cm,
Weight: 95 kg,
Valid inspection and ready to use,
Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international market,
Financing options (Poland only): Installments, Leasing, Loan,

