



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



FRESENIUS VOLUMAT AGILIA PL infusion pump (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in France,
Power supply: 230 V,
Backup power supply: via 2 new built-in batteries,
Frequency: 50/60 Hz,
Power: 15 VA,
Technical specifications:
Flow rate range:
Normal mode:
1-1500 ml/h (in 1 ml/h increments),
Micro mode:
0.1-99.9 ml/h (in 0.1 ml/h increments),
Flow rate can be limited depending on the drug name (soft and hard limits) using Vigilant®, a solution ensuring the safety of intravenous medications,
Flow rate accuracy: +/- 5%
Infusion modes:
ml/h modes: volume + flow rate,
Volume + Time,
Flow Rate + Time,
Volume + Time + Rate,
Rate Only (if optional drip sensor is used),
Ramp/Drop,
Sequential/Intermittent,
Secondary/Piggyback,
Tramp/Min,
Dose Modes:
ng/h, ng/kg/min, ng/kg/h, ng/m²/24h, µg/min, µg/h, µg/kg/min, µg/kg/h, µg/m²/min, µg/m²/h, mg/min, mg/h, mg/24h, mg/kg/min, mg/kg/h, mg/m²/24h, mg/m²/min, mg/m²/h, g/h, g/kg/min, g/kg/h, g/m²/min, g/m²/24h mU/min, mU/h, mU/kg/min, mU/kg/h, mU/m²/min, mU/m²/h, U/min, U/h, U/kg/h, U/kg/24h, U/m²/min, U/m²/h, mEq/min, Eq/h, mEq/kg/min, mEq/kg/h, mEq/m²/min, mEq/m²/h, mEq/m²/24h, mmol/h, kcal/h, kcal/24h, kcal/kg/h, kcal/kg/24h. With or without a loading dose.
Dose Time Mode:
Dose calculation is based on the dose and dilution over a predetermined time.
OCS Patent:
The Occlusion Check System is the only self-test that verifies the proper operation of the pump and its kit, thus preventing the risk of free flow.
Bolus: 1500 ml/h, adjustable from 50 ml/h to 1500 ml/h in 50 ml/h increments. Manual or programmable,
Programming modes:
No drug names: infusion without displaying drug names,
Drug labeling: infusion with displayed drug names,
Vigilant® Drug'Lib: IV administration mode secured with soft and hard limits,
Customized drug library with Vigilant, including an IV drug safety solution,
2 pressure modes: variable (from 50 to 750 mmHg, in 50 mmHg increments) or 3 preset levels,
DPS: Dynamic Pressure System (DPS) alerts you to pressure changes. Therefore, the risk of blockage or possible leakage in the infusion line can be predicted.
Menu language selection: FR, ENG, GER, NR, SE, DU, GR, IT, ESP, PL, TUR, RUS, RUM, CS, HUN, DUN.
Includes:
Power cord,
Mounted holder for hanging on a stand,
Instruction manual in Polish.
Dimensions: 13 x 13.5 x 14.5 cm.
Weight: 2 kg.
Updated inspection and ready for operation.
Issued Technical Passport (Service Report) valid for 12 months.
Warranty:
6 months for the domestic market (Poland),
3 months for the international market.
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international market.
Financing options (Poland only): Installments, Leasing. Loan,

