



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MGE WBRKT ambulance holder for ZOLL R-SERIES DEFIBRILLATOR (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- UK Production, (2013),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Ambulance handle designed for **the ZOLL R-SERIES** defibrillator ,
- It has a built-in 230 V power connector,
- Solid and durable construction ensures safety and comfort of use,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

