



DRAGER MOVA CART endoscopic trolley (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Made in Germany,

Specifications:

Sliding drawer dimensions: 50 x 40 cm,

Double-sided handles for easy transport,

Mobile trolley on 4 wheels (2 with brakes),

Dimensions:

Shelves: 45 x 60 cm,

Sliding drawer dimensions: 50 x 40 cm,

Total dimensions: 80 x 55 x 142 cm,

Maximum shelf load capacity: 50 kg,

Weight: 35 kg,

Updated inspection and ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility to extend the warranty for up to 12 months of subsidy for domestic and international markets. Financing options (Poland only):

Installments, Leasing, Loan.

